




Care & Cure Surgico

"Our Creation Is For Your Health"



LARYNGOSCOPE
ENT DEVICES
CATALOGUE


www.carecuresurgico.com

Company Profile

Care & Cure Surgico is dedicated to provide you the personalized service you deserve. We have been supplying the medical industry with the finest in Surgical Instruments , Dental instruments , Laryngoscope, ENT Devices , Hospital Furniture and Stainless Steel Hollowware for all specialties for over 40 years. Our office is located in Sialkot, Pakistan, we offer complete in-house repair service and custom manufacturing in addition to our already broad line of products. Care & Cure Surgico is a one stop shop for many clients it also assist turn key projects, works with MOH and NGOs, Hospitals and present around the world with large number of the distributors worldwide

When you buy from Care & Cure Surgico, you will enjoy the services and quality we aim to deliver from years. Care & Cure Surgico comes with not only the best products, but the best services as well. If you ever have a question, We're here to help. We will be more than happy to give you the support you need.

Certifications & Standards Details



To Our Existing Customers

We would like to thank you for continuing to stay with the best. We appreciate your business.

To Prospective Customer

We look forward to starting a lasting relationship with you.

Guarantee

Care & Cure Surgico products are guaranteed to be free of defects in materials and workmanship. Care & Cure Surgico products comes with guarantee or warranty as per the products category but usually as per international guarantee rules. Any product that is determined to be defective will be either repaired or replaced, at our discretion. Product not properly cared for or used for purposes not intended will not be covered. Care & Cure Surgico liability under this guarantee is limited to repair or replacement of defective items.

OEM Services

We have been providing OEM services to various bands in the health care industry all over the world for the past 25 years

OEM SERVICES PROVIDED, WITH THE UTMOST DISCRETION

Care & Cure Surgico

"Our Creation Is For Your Health"

www.carecuresurgico.com



www.carecuresurgico.com

Care & Cure Surgico

"Our Creation Is For Your Health"

Mc – Intosh Laryngoscopes

Standard Illumination

Adult Size



20-200-15

Handle Standard Illumination

Size: Adult

Handle options: Child & Stubby

Standard Sets

Contents	Code
■ Set of 3 Blades Size 0, 1, 2 and Adult Handle	■ 20-200-30
■ Set of 3 Blades Size 1, 2, 3 and Adult Handle	■ 20-200-31
■ Set of 4 Blades Size 0, 1, 2, 3 and Adult Handle	■ 20-200-40
■ Set of 4 Blades Size 1, 2, 3, 4 and Adult Handle	■ 20-200-41

Sets available in hard case or soft pouch



20-100-00

Mc - Intosh Blade Standard Illumination


Size: 0



20-100-01

Mc - Intosh Blade Standard Illumination

Size: 1



20-100-02

Mc - Intosh Blade Standard Illumination

Size: 2



20-100-03

Mc - Intosh Blade Standard Illumination

Size: 3




20-100-04

Mc - Intosh Blade Standard Illumination

Size: 4

Specification:

- Manufactured From AISI 304 Stainless Steel
- All Blades and Handles Compatible with ISO 7376 Standard
- Autoclavable
- Available in Halogen and LED Illumination
- Bulbs 2.5V
- Handles requires 2 x C size batteries
- Hard Case and Soft Pouch available



Head



Handle



Blade



Mc – Intosh Laryngoscopes

Fiber Optic Illumination

Adult Size



21-200-15

Handle Fiber Optic Illumination
Size: Adult

Handle options: Child & Stubby

Standard Sets


Contents	Code
▪ Set of 3 Blades Size 0, 1, 2 and Adult Handle	▪ 21-200-30
▪ Set of 3 Blades Size 1, 2, 3 and Adult Handle	▪ 21-200-31
▪ Set of 4 Blades Size 0, 1, 2, 3 and Adult Handle	▪ 21-200-40
▪ Set of 4 Blades Size 1, 2, 3, 4 and Adult Handle	▪ 21-200-41

Sets available in hard case or soft pouch



21-100-00

Mc - Intosh Blade Fiber Optic Illumination
Size: 0



21-100-01

Mc - Intosh Blade Fiber Optic Illumination
Size: 1



21-100-02

Mc - Intosh Blade Fiber Optic Illumination
Size: 2



21-100-03

Mc - Intosh Blade Fiber Optic Illumination
Size: 3




21-100-04

Mc - Intosh Blade Fiber Optic Illumination
Size: 4

Specification:

- Fiber Optic Handle with 3.5mm dia for improved light transmission and longer life.
- Manufactured From AISI 304 Stainless Steel
- All Blades and Handles Compatible with ISO 7376 Standard
- Autoclavable
- Available in Halogen and LED Illumination
- Bulbs 2.5V
- Handles requires 2 x C size batteries
- Hard Case and Soft Pouch available



Head



Handle



Blade



www.carecuresurgico.com

Care & Cure Surgico

"Our Creation Is For Your Health"

Miller Laryngoscopes Standard Illumination Child / Infant Size



22-200-15

Handle Standard Illumination

Size: Child / Infant

Handle options: Adult & Stubby

Standard Sets

Contents	Code
■ Set of 3 Blades Size 00, 0, 1 and Child Handle	■ 22-200-30
■ Set of 3 Blades Size 0, 1, 2 and Child Handle	■ 22-200-31
■ Set of 4 Blades Size 00, 0, 1, 2 and Child Handle	■ 22-200-40
■ Set of 4 Blades Size 0, 1, 2, 3 and Child Handle	■ 22-200-41

Sets available in hard case or soft pouch



22-100-00

Miller Blade Standard Illumination

Size: 00



22-100-01

Miller Blade Standard Illumination

Size: 0



22-100-02

Miller Blade Standard Illumination

Size: 1



22-100-03

Miller Blade Standard Illumination

Size: 2



22-100-04

Miller Blade Standard Illumination

Size: 3

Specification:

- Manufactured From AISI 304 Stainless Steel
- All Blades and Handles Compatible with ISO 7376 Standard
- Autoclavable
- Available in Halogen and LED Illumination
- Bulbs 2.5V
- Handles requires 2 x AA size batteries
- Hard Case and Soft Pouch available



Head



Handle



Blade

Miller Laryngoscopes

Fiber Optic Illumination

Adult Size

www.carecuresurgico.com

Care & Cure Surgico

"Our Creation Is For Your Health"



23-200-15

Handle Fiber Optic Illumination

Size: Adult

Handle options: Child & Stubby

Standard Sets

Contents	Code
■ Set of 3 Blades Size 0, 1, 2 and Adult Handle	■ 23-200-30
■ Set of 3 Blades Size 1, 2, 3 and Adult Handle	■ 23-200-31
■ Set of 4 Blades Size 0, 1, 2, 3 and Adult Handle	■ 23-200-40
■ Set of 4 Blades Size 1, 2, 3, 4 and Adult Handle	■ 23-200-41

Sets available in hard case or soft pouch



23-100-00

Miller Blade Fiber Optic Illumination

Size: 0



23-100-01

Miller Blade Fiber Optic Illumination

Size: 1



23-100-02

Miller Blade Fiber Optic Illumination

Size: 2



23-100-03

Miller Blade Fiber Optic Illumination

Size: 3




23-100-04

Miller Blade Fiber Optic Illumination

Size: 4

Specification:


- Fiber Optic Handle with 3.5mm dia for improved light transmission and longer life.
- Manufactured From AISI 304 Stainless Steel
- All Blades and Handles Compatible with ISO 7376 Standard
- Autoclavable
- Available in Halogen and LED Illumination
- Bulbs 2.5V
- Handles requires 2 x C size batteries
- Hard Case and Soft Pouch available



Head



Handle



Blade



www.carecuresurgico.com

Care & Cure Surgico

"Our Creation Is For Your Health"

Mc – Intosh Integrated Laryngoscopes

Fiber Optic Illumination

Adult Size



24-200-15

Handle Fiber Optic Illumination
Size: Adult

Handle options: Child & Stubby

Standard Sets

Contents	Code
▪ Set of 3 Blades Size 00, 0, 1 and Adult Handle	▪ 24-200-30
▪ Set of 3 Blades Size 1, 2, 3 and Adult Handle	▪ 24-200-31
▪ Set of 4 Blades Size 0, 1, 2, 3 and Adult Handle	▪ 24-200-40
▪ Set of 4 Blades Size 1, 2, 3, 4 and Adult Handle	▪ 24-200-41

Sets available in hard case or soft pouch



24-100-00

Mc - Intosh Blade Fiber Optic Illumination
Size: 00



24-100-01

Mc - Intosh Blade Fiber Optic Illumination
Size: 0



24-100-02

Mc - Intosh Blade Fiber Optic Illumination
Size: 1



24-100-03

Mc - Intosh Blade Fiber Optic Illumination
Size: 2



24-100-03

Mc - Intosh Blade Fiber Optic Illumination
Size: 3




24-100-04

Mc - Intosh Blade Fiber Optic Illumination
Size: 4

Specification:

- Fiber Optic Handle with 3.5mm dia for improved light transmission and longer life.
- Manufactured From AISI 304 Stainless Steel
- All Blades and Handles Compatible with ISO 7376 Standard
- Autoclavable
- Available in Halogen and LED Illumination
- Bulbs 2.5V
- Handles requires 2 x C size batteries
- Hard Case and Soft Pouch available



Head



Handle



Blade



McCoy Flexi - Tip Laryngoscopes

Fiber Optic Illumination

Adult Size



25-200-15
Handle Fiber Optic Illumination
Size: Adult



25-100-00
McCoy Blade Fiber Optic Illumination
Size: 1



25-100-02
McCoy Blade Fiber Optic Illumination
Size: 2



25-100-02
McCoy Blade Fiber Optic Illumination
Size: 3




Standard Sets

Contents	Code
■ Set of 2 Blades Size 2, 3 and Adult Handle	■ 25-200-30
■ Set of 2 Blades Size 3, 4 and Adult Handle	■ 25-200-31
■ Set of 3 Blades Size 2, 3, 4 and Adult Handle	■ 25-200-40

Sets available in hard case or soft pouch

Specification:

- Flexible Tip for difficult Intubation
- Fiber Optic Handle with 3.5mm dia for improved light transmission and longer life.
- Manufactured From AISI 304 Stainless Steel
- All Blades and Handles Compatible with ISO 7376 Standard
- Autoclavable
- Available in Halogen and LED Illumination
- Bulbs 2.5V
- Handles requires 2 x C size batteries
- Hard Case and Soft Pouch available



Head



Handle



Blade



www.carecuresurgico.com

Care & Cure Surgico

"Our Creation Is For Your Health"

McCoy Flexi - Tip Integrated Laryngoscopes

Fiber Optic Illumination

Adult Size



26-200-15

Handle Fiber Optic Illumination
Size: Adult



26-100-00

McCoy Blade Fiber Optic Illumination
Size: 1



26-100-02

McCoy Blade Fiber Optic Illumination
Size: 2



26-100-02

McCoy Blade Fiber Optic Illumination
Size: 3




Standard Sets

Contents	Code
■ Set of 2 Blades Size 2, 3 and Adult Handle	■ 26-200-30
■ Set of 2 Blades Size 3, 4 and Adult Handle	■ 26-200-31
■ Set of 3 Blades Size 2, 3, 4 and Adult Handle	■ 26-200-40

Sets available in hard case or soft pouch

Specification:

- Flexible Tip for difficult Intubation
- Fiber Optic Handle with 3.5mm dia for improved light transmission and longer life.
- Manufactured From AISI 304 Stainless Steel
- All Blades and Handles Compatible with ISO 7376 Standard
- Autoclavable
- Available in Halogen and LED Illumination
- Bulbs 2.5V
- Handles requires 2 x C size batteries
- Hard Case and Soft Pouch available



Head



Handle



Blade

ENT Devices Fiber Optic Illumination

www.carecuresurgico.com

Care & Cure Surgico
"Our Creation Is For Your Health"



27-200-15

Brass Otoscope Fiber Optic
Brass Otoscope Conventional



28-200-15

ENT Set Complete



29-200-15

Brass Ophthalmoscope

Specification:

- Manufactured From Brass Chrome plated
- 3x Magnifying glass
- Available in Halogen and LED Illumination
- Bulbs 2.5V
- Handles requires 2 x C size batteries
- Hard Case and Soft Pouch available
- Fiber Optic Light Source



www.carecuresurgico.com

Care & Cure Surgico

"Our Creation Is For Your Health"

Mini Otoscope Fiber Optic Illumination



30-200-15

Mini Otoscope Fiber Optic
Complete Set of Disposable Cannulas



Color Options

Specification:

- Manufactured From Plastic Material
- 3x Magnifying glass
- Available in Halogen and LED Illumination
- Bulbs 2.5V
- Handles requires 2 x AA size batteries
- Soft Pouch available in matching colors
- Fiber Optic Light Source



Standard Accessories

Care & Cure Surgico

"Our Creation Is For Your Health"

Nasal Speculum



31-200-15



32-200-15



33-200-15

Bulbs



34-200-15



35-200-15



36-200-15



37-200-15

Ear Speculas



38-200-15



39-200-15



40-200-15



41-200-15

Ear Speculas & Testing Tube



42-200-15



43-200-15



44-200-15



www.carecuresurgico.com

Care & Cure Surgico

"Our Creation Is For Your Health"

Packaging Options



45-200-15



46-200-15



47-200-15



48-200-15

Please, flash this code of
Our Catalogue



If required, we will gladly send you the Catalogue of our others products.



Care & Cure Surgico P.O Box 919, Noal More,
Roras Road, Sialkot - 51310 Pakistan.

T: +92 52 3563200 T: +92 52 3252500

F: +92 52 3553777

E: info@carecuresurgico.com

W: www.carecuresurgico.com

export.ccs



www.facebook.com/carecuresurgico

www.twitter.com/carecuresurgico

www.linkedin.com/company/care-&-cure-surgico

www.plus.google.com/+Carecuresurgico



ISO 9001: 2008
ISO 13485: 2003

