



Care & Cure Surgico

"Our Creation Is For Your Health"



Dental Cassettes

"Scan these QR codes to Reach us"



<https://www.linkedin.com/company/care-&-cure-surgico>



<https://www.twitter.com/carecuresurgico>



https://www.instagram.com/care_and_cure_surgico/



<https://www.facebook.com/carecuresurgico>



+923458420345





Care & Cure Surgico

"Our Creation Is For Your Health"

Company Overview:

Care & Cure Surgico is dedicated to provide you the personalized service you deserve. We have been supplying the medical industry with the finest in Surgical Instruments , Dental instruments , Laryngoscope, ENT Devices, Hospital Furniture and Stainless Steel Hollowware for all specialties for over 40 years. Our office is located in Sialkot. Pakistan, we offer complete in-house repair service and custom manufacturing in addition to our already broad line of products. Care & Cure Surgico is a one stop shop for many clients it also assist turn key projects, works with MOH and NGOs, Hospitals and present around the world with large number of the distributors worldwide When you buy from Care & Cure Surgico, you will enjoy the services and quality we aim to deliver from years. Care & Cure Surgico comes with not only the best products, but the best services as well. If you ever have a question, We're here to help. We will be more than happy to give you the support you need.

Our Guarantee:

Care & Cure Surgico products are guaranteed to be free of defects in materials and workmanship. Care & Cure Surgico products comes with guarantee or warranty as per the products category but usually as per international guarantee rules. Any product that is determined to be defective will be either repaired or replaced, at our discretion. Product not properly cared for or used for purposes not intended will not be covered. Care & Cure Surgico liability under this guarantee is limited to repair or replacement of defective items.

Contact us:



Roras Road Sialkot - 51310 Pakistan



+(92) 523563200



info@carecuresurgico.com



<https://carecuresurgico.com/>



Care & Cure Surgico

"Our Creation Is For Your Health"

Design & Purpose

Our Products are designed to consistently meet the high demand of our customers.

we believe in quality and service and we aren't going to compromise when it comes to meeting the justified Expectations of our customers.

Our aim is to provide quality product for Health care sector around the globe to achieve the goal of healthy human being

R&D:


The Passion and desire to make a surgeries as smooth and as easy our department continues to develop smooth products that matter in health sector today

This will of Continuous Struggles Comes from our relationship as a human being. We are committed to bring out the best solution for our customers

OEM Services

We have been providing OEM/Bulk services to various bands in the health care industry all over the world for the past 25 years

Certificates:



To Our Existing Customers

We would like to thank you for continuing to stay with the best. We appreciate your business.

To Prospective Customer

We look forward to starting a lasting relationship with you.

"Scan these QR codes to Reach us"



<https://www.linkedin.com/company/care-&-cure-surgico>



<https://www.twitter.com/carecuresurgico>



+923458420345



https://www.instagram.com/care_and_cure_surgico/



<https://www.facebook.com/carecuresurgico>

Copyright law



All rights reserved. The contents of this catalog, including but not limited to text, images, and product descriptions, are protected by copyright laws and belong to Care & Cure Surgico. Reproduction, distribution, or unauthorized use of any materials from this catalog, in whole or in part, is strictly prohibited without the express written consent of Care & Cure Surgico.

Care & Cure Surgico retains the exclusive rights to the design, layout, and compilation of this catalog. The trademarks, logos, and brand names featured in this catalog are the property of Care & Cure Surgico and are protected under national and international trademark laws. Any unauthorized use, reproduction, or imitation of these trademarks is strictly prohibited.

Care & Cure Surgico takes intellectual property rights seriously and will take legal action against any infringement or misuse of its copyrighted materials, trademarks, or intellectual property.

For inquiries regarding the use of materials or obtaining permission to reproduce any part of this catalog, please contact Care & Cure Surgico's legal department.

This catalog serves as a comprehensive reference for your professional needs. Please note that not all items available in our inventory are included in this catalog. Care & Cure Surgico is dedicated to providing you with complete technical and commercial information as required.

Care & Cure Surgico operates under a Certified Quality Management System, adhering to ISO 9001:2015 Quality Standards and UNI EN ISO 13485:2016. Our commitment ensures that all medical devices featured in this catalog fully comply with current EEC directives governing medical devices.


© 1973-2023 Care & Cure Surgico , All Rights Reserved , Copyright Ordinance 1962



INSTRUMENT CASSETTES

Organized Procedure Kits save time and improve safety! We are pleased to introduce a new line of cassette trays to organize and protect instruments during handling and transportation. Removable, silicone instrument holders afford you the option of an open configuration in any cassette. Cassettes with detachable lids feature moveable instrument holders to customize instrument configuration. Choose from various cassette sizes and styles to fit your needs.

- High quality, heavy gauge stainless steel.
- Optimum slot position for effective cleaning and drainage.
- Easy open/closing latch system.
- Silicone Inserts for secure instrument placement.



SMALL CASSETTES

Typical Applications Include:

- Diagnostic Kit
- Root Planning Kit
- New Patient Kit
- Composite Kit

PRODUCT NUMBER	SIZE	CONFIGURATION	FEATURES
CAS-100-1	7" x 2¾" x 11/3"	5 Instruments	detachable Lid/Instrument Holders
CAS-100-2	7" x 21/3" x 6/7"	Instruments	Fixed Lid/Instrument Holders



MEDIUM CASSETTES

Typical Applications Include:

- Hygiene Kit
- Restoration Kit
- Cavity Prep Kit
- Endo Kit
- Crown & Bridge Kit

PRODUCT NUMBER	SIZE	CONFIGURATION	FEATURES
CAS-100-3	7"x41/3"x11/3"	7Instruments	detachable Lid/Instrument Holders
CAS-100-4	7" x 3½" x 6/7"	7 Instruments	Fixed Lid/Instrument Holders



MEDIUM CASSETTES, CONTINUED

PRODUCT NUMBER	SIZE	CONFIGURATION	FEATURES
CAS-100-5	7"x5¼"x11/3"	10 Instruments	detachable Lid/Instrument Holders
CAS-100-6	7" x 5" x 6/7"	10 Instruments	Fixed Lid/Instrument Holders



LARGE CASSETTES

Typical Applications Include:

- Periodontal Surgery
- Restoration Kit
- Extraction

PRODUCT NUMBER	SIZE	CONFIGURATION	FEATURES
CAS-100-7	11½"x7"x11/3"	20 Instruments	detachable Lid/Instrument Holders
CAS-100-8	11" x 7¼" x 1"	20 Instruments	Fixed Lid/Instrument Holders



EXTRA LARGE CASSETTE

Available with specially designed silicone inserts to hold up to 12 extracting forceps. or choose the open configuration recommended for large surgical instrument sets.

PRODUCT NUMBER	SIZE	CONFIGURATION	FEATURES
CAS-100-9	15" x 8¾" x 1½"	12 Forceps	detachable Lid/Instrument Holders
CAS-100-10	15" x 8¾" x 1½"	open	detachable Lid






INSTRUMENT HOLDERS


univerSAl

PRODUCT NUMBER	SIZE	CONFIGURATION	FEATURES
CAS-100-11	5"	10 Instruments	3 pc Silicone Set, blue
CAS-100-12	10½"	20 Instruments	3 pc Silicone Set, blue



CLIP IN CASSETTES


PRODUCT NUMBER	SIZE	CONFIGURATION	FEATURES
CAS-100-13	7" x 2⅔" x 6/7"	5 Instruments	open top and bottom



PRODUCT NUMBER	SIZE	CONFIGURATION	FEATURES
CAS-100-14	7" x 3½" x 1"	10 Instruments	open top with Instrument Holders



Material we use in these Sets



Stainless Steel

Surgical Dental tools must be made of high-quality materials to ensure the highest levels of performance and safety.

Dental Cassettes need to be durable, hygienic, and, depending on the surgery, lightweight and unobtrusive.

Needle Holders are the norm. These Dental Cassettes offer excellent strength and can be sterilized and reused many times.

Medical- or surgical-grade steels are typically stainless steels composed of chromium, carbon, nickel, and molybdenum.


The property of being alloys makes them more rigid and more corrosion-resistant.

Common examples include austenitic 316 and martensitic 420 and 440. Conventionally, these are widely used materials for metal retractors used in surgery.

Stainless Steel: Stainless steel cassettes are durable, resistant to corrosion, and can withstand the rigors of sterilization processes. They are often preferred for their strength and longevity. Stainless steel cassettes are commonly used in dental practices that require frequent sterilization and heavy instrument usage.

They offer excellent protection for instruments and are easy to clean and maintain.


Quality Goal





What we do after customer send us inquiry

SUMARIZING






Care & Cure Surgico

"Our Creation Is For Your Health"

WHERE WE LOCATED?

Care & Cure Surgico, Roras Road



Care & Cure Surgico
4.4 ★★★★★ 8 reviews
Medical equipment manufacturer


Directions Save Nearby Send to phone Share

Noal More, Roras Rd, Sialkot, Punjab 51310, Pakistan

سیالکوٹ، ضلع سیالکوٹ، پنجاب، Noal More, Roras Rd 51310

Open · Closes 5 PM

carecuresurgico.com



Noal More, Roras Road Sialkot - 51310 Pakistan



Care & Cure Surgico

"Our Creation Is For Your Health"


THANK YOU

We look forward to working with you

CONTACT US:

 Roras Road Sialkot - 51310

Pakistan

 +(92) 523563200

 info@carecuresurgico.com

"Scan these QR codes to Reach us"



<https://www.twitter.com/carecuresurgico>



<https://www.facebook.com/carecuresurgico>



<https://www.linkedin.com/company/care-&-cure-surgico>



https://www.instagram.com/care_and_cure_surgico/



"Please Scan QR CODE TO CHAT WITH US ON WHATSAPP"



+923458420345

